

Information is specific to use of each pump's AID features

CALCULATE

Auto correction boluses (max. every

Auto corrections can be turned on or

Auto correction boluses (max once/hr)

No automated boluses. Algorithm will

increase basal doses up to 400% to

No automated boluses. Algorithm will

increase basal doses up to the max

basal rate programmed in the pump

All meal bolus doses and correction

minutes as needed if glucose > CGM

ADJUST

• Can users

adjust the

correction

bolus

target?

No, fixed at

120 mg/dL

No, fixed at

110 mg/dL

Yes, same

Glucose")

Yes, same

target

Range")

as algorithm

("Correction

Yes, same

as algorithm

target ("CGM

Target")

REVERT

EDUCATE

already decreased or suspended insulin delivery and treating with too many carbs may result in large rebound hyperglycemia.

Consider treating mild hypoglycemia with less carbohydrates (5-10 g) than the traditional rule of 15g. If hypoglycemia occurs, the algorithm will have

Consider using the 100 mg/dL Target and Active Insulin Time of 2 hours for optimal system performance as long as hypoglycemia is not >4%.

The sensor glucose value auto-populates into the bolus menu for correction bolus calculation. SmartGuard will adjust the bolus dose based on the CGM

Do not enter "fake carbs" to try to get more insulin from the system. This will result in an increased risk of hypoglycemia, and greater glucose variability.

It is best NOT to override the bolus calculator's suggested dose (although there may be exceptions). The bolus calculator will subtract IOB from increased

Control-IQ allows programming of more than 1 personal profile, where different basal rates, carb ratios and correction factors can be used. Additional profile

Tap "Use Sensor" to add the sensor glucose value and trend into the bolus. The bolus calculator may adjust the recommended correction bolus dose based

It is best NOT to override the bolus calculator's suggested dose (although there may be exceptions). The bolus calculator will subtract IOB from increased

Insulin suspension may occur if glucose is trending down, even if the glucose level is above the programmed Target Glucose. This is expected and will be

Use the Pre-Meal Preset up to 1 hour before meals to help reduce post-prandial glucose spikes. Pre-Meal Preset is a lower correction range allowing

Indicate expected carbohydrate absorption for meals to better match insulin delivery to food type and help algorithm track active carbohydrates more

User should "announce" meals at the start of the meal by indicating meal type ("Breakfast" "Lunch" or "Dinner") and meal size ("Usual for me", "More" than

The iLet is designed to automate all insulin delivery, and continuously adapts with no user interaction except for meal announcements. Users cannot give

Don't use meal announcements to try to correct high glucose levels; this will disrupt the system's adaptation and increase the chance of hypoglycemia.

It's important to carry a BG meter at all times so the user can use the BG-run mode if there are unexpected problems with the CGM at any time. BG-run mode

data)

Can others see data remotely?

CareLink Connect app (pump + CGM

Dexcom Follow app (CGM data)

Dexcom Follow app (CGM data)

*If using Freestyle Libre 2 Plus, there

is no option for remote data sharing

twiist insight app for remote pump

Dexcom Follow app (CGM data)

Bionic Circle app (pump + CGM data)

and sensor data sharing

*If using Freestyle Libre 2 Plus, there

is no option for remote data sharing

Is data automatically stored in

Automatic uploads to CareLink via

Automatic uploads to Source via t:connect mobile app or via Mobi app

Automatic uploads to Glooko or

Automatic uploads to Tidepool Data

Automatic uploads to the Beta Bionics

Platform

portal via iLet app

Discover after linking device

the cloud?

MiniMed mobile app

usual, or "Less" than usual) relative to the user's typical carbohydrate intake for each meal type. Do not announce the meal if > 30 min. after eating.

In the first week of using iLet, space meals at least 4 hours apart and eat primarily "Usual for me" meals to help the iLet learn meal bolus doses.

SENSOR/SHARE

Can user see real-time data on

MiniMed mobile app (pump + CGM

personal cell phone?

Dexcom G6/G7 mobile app

t:slim = t:Connect mobile app

Mobi = Mobi app (pump + CGM data)

Omnipod 5 app (pump + CGM data,

availability of app varies by region)

also used to operate pump;

Dexcom G6/G7 mobile app

twiist app (pump + CGM data)

Dexcom G6/G7 mobile app

iLet Mobile app (pump + CGM data)

(CGM data)

(pump + CGM data)

• Emoji options: Lollipop: Fast (30 min), Taco: Medium (3 hours), Pizza: Slow (5 hours); or manually enter absorption time 30 min-8 hours

Pump button available to bolus on pump. One press = 1 unit. Can give max of 10 unit bolus via pump button.

lasts max. 72 hours, so if CGM wear is interrupted for longer, a back up insulin delivery plan is vital.

data)

(CGM data)

(CGM data)

Consider turning the reverse correction OFF. The reverse correction will reduce the meal bolus dose if the glucose level is below the target glucose.

as algorithm

target ("Target

bolus doses are automated. Auto

correction boluses max. every 5

Can users

adjust

active

time?

insulin

No, fixed at

5 hrs when

Yes

No, fixed

at 6 hours

Control-IQ is

target setting.

Can users

correction

(sensitivity)?

adiust

factor

No (the

programmed

sensitivity is

not used for

calculations

when in SmartGuard)

Yes

Yes

Yes

There are no pump settings programmed into the iLet.

continuously adapted without the use of any programmed

Is there a limited automation mode the system may revert to if

there is a loss of CGM communication or other reasons?

Yes, Safe Basal: the pump will deliver a basal rate determined by the

User needs to enter a BG value into the pump before the "time to exit"

algorithm, but without glucose-dependent basal adjustments and no auto

May activate due to max/min insulin delivery constraints, loss of CGM data

If there is loss of CGM data, the pump will deliver the programmed basal rates

without glucose-dependent basal adjustments and no auto correction boluses

Yes, Automated Limited: the Pod will deliver a basal rate determined by the

1. If no CGM data for ≥ 20 min. Pod will resume full insulin automation once

2. If there is an "Automated Delivery Restriction" alarm (if insulin has been suspended too long or if max delivery too long). Will remain in Automated

Yes, BG-run mode: If the iLet loses communication with the CGM, it will

prompt the user to enter BG values periodically. As long as the user enters

BG values into the iLet, it will continue to automate all insulin delivery based

on the entered BG values and previously stored information on the user's

The user can continue to announce meals in BG-run mode to receive meal

Pre-bolus for all meals and snacks, ideally 10-15 min before eating.

Follow system prompts for "BG Required" to stay in SmartGuard.

Pre-bolus for all meals and snacks, ideally 10-15 min before eating.

Pre-bolus for all meals and snacks, ideally 10-15 min before eating.

more aggressive basal insulin leading up to mealtime only.

Apple Watch compatibility – can bolus via Apple watch.

a manual bolus. A hands-off approach is necessary.

Which CGM is compatible?

Guardian 4

Simplera Sync

*CGM options may vary by region

Dexcom G6 and G7: Use of Dexcom

cannot use the Dexcom receiver when

G6 or G7 mobile app is optional;

the Dexcom is paired to the pump.

Freestyle Libre 2 Plus (t:slim only):

Must connect CGM to the pump via t:connect Mobile app. Cannot use Freestyle Libre apps or reader.

Dexcom G6 and G7: Must use Dexcom

Cannot use the Dexcom receiver when

Freestyle Libre 3+: Must pair on twiist

app, cannot use Libre apps or reader

Dexcom G6 and G7: Use of Dexcom

cannot use the Dexcom receiver when

G6 or G7 mobile app is optional;

the Dexcom is paired to the iLet.

Libre apps or reader.

Freestyle Libre 3 Plus: Must pair on iLet app only; cannot use Freestyle

G6/G7 mobile app (on personal cell

phone) to use Automated Mode.

the Dexcom is paired to the Pod.

Freestyle Libre 2 Plus: Must start sensor on Omnipod 5 controller; cannot use Freestyle Libre apps or

reader.

automated insulin delivery, helping to reduce the chance of hypoglycemia.

Mobi: requires Mobi app on personal cell phone (iOs only) for pump control.

automated insulin delivery, helping to reduce the chance of hypoglycemia.

short in duration (e.g. 5-15 min) if the glucose level does not continue to drop. Wear Pod and Dexcom in "line of sight" to optimize Bluetooth communication.

twiist app only compatible with iOs. Must have an iPhone to use twiist system.

value and insulin on board. The user is not able to change or override the suggested dose.

The sensor glucose value auto-populates into the bolus menu for correction bolus calculation.

t:slim X2: can give bolus doses remotely from a cell phone when using the t:connect mobile app.

to help with changing insulin needs (e.g., menstrual cycle, illness, long sporting events, etc.).

Program the sleep schedule to ensure Sleep Activity activates each day automatically.

iLet can operate in BG-run mode for up to 72 hrs.

No. If there is loss of CGM data for more than 15 minutes, the pump will revert

algorithm, but without glucose-dependent basal adjustments. May activate

All insulin delivery is automated by the algorithm and

or system concerns about sensor accuracy.

expires to prevent SmartGuard exit.

No, there is no limited automation mode.

Limited until the user clears the alarm.

to the programmed basal rates.

correction bolus

to help correct hyperglycemia.

help correct hyperglycemia.

if glucose is predicted to be >180

mg/dL in 30 min.

5 min) if glucose is >120 mg/dL.

Algorithm target glucose/target

3 target options: 100, 110, 120 mg/dL

Target range: 112.5-160 mg/dL

5 target options: 110, 120, 130, 140,

Can set multiple target settings

Target range (called "Correction

87-180 mg/dL. Algorithm targets

Range"), can set any range between

Can set multiple ranges throughout

3 target options: Usual (120 mg/dL),

User can set up to 2 target settings

What are the special features

in automated insulin delivery?

Temp Target: Changes target glucose

delivery for chosen duration (30 min -

to 150 mg/dL to reduce auto-basal

24 hr) and disables auto correction

Exercise Activity: Changes target

range to 140-160 mg/dL to reduce

Sleep Activity: Narrows target range to 112.5-120 mg/dL and disables auto

Can program a sleep schedule or

Activity Feature: Changes target

basal delivery for chosen duration

Pre-Meal Preset (correction range

for up to 1 hour before meals, 67-110

mg/dL, typically set lower than general

Workout Preset (set a correction range

retrospectively (algorithm tracks active carbohydrates to use for glucose predictions. Can correct or change carbohydrate entries as needed).

to use for exercise, 87-250 mg/dL).

Pause Insulin Feature: Users can

a specified timeframe.

mode.

pause (suspend) insulin delivery for

no insulin dose automation)?

an exit to manual mode.

Modify Carbohydrate entries

glucose to 150 mg/dL and decreases

the doses by ~50% to reduce adaptive

Lower (110 mg/dL), Higher (130

throughout 24 hr period

middle of the range.

24 hour day.

mg/dL)

boluses.

basal delivery.

correction boluses.

manually start/stop.

(1-24 hrs).

correction range).

per 24 hr period

range?

150 mg/dL

Which insulin does the user

User gives boluses for meals by

bolus menu / bolus calculator.

entering total grams of carbs in the

User can deliver correction boluses

User gives boluses for meals by

bolus menu / bolus calculator.

entering total grams of carbs in the

User can deliver correction boluses

User gives boluses for meals by

bolus menu / bolus calculator.

entering total grams of carbs in the

User can deliver correction boluses

User gives boluses for meals by

entering total grams of carbs in the

User can deliver correction boluses

"announcement" to prompt the iLet to

deliver a meal bolus, which involves

estimating the carbohydrate amount

for each meal ("Usual for Me"/"More"

Which pump settings can be

insulin delivery (automated

Auto Basal Target: 100, 110, 120

Active Insulin Time (2 hrs for most

5 target options: 110, 120, 130, 140,

Maximum Basal Rate (used as

Correction Range. Consider a

max basal for basal automation).

Recommended to set at 3-4x highest

20mg/dL range (e.g., 90-110 mg/dL).

User can set up to 2 target settings

per 24 hr period: "CGM Target" and

"Secondary CGM Target".

mg/dL; only 1 target can be set.

insulin delivered by the

aggressive insulin delivery)

algorithm)?

Basal rates

Correction factor

Target Glucose

150 mg/dL

Basal rates

basal rate.

When will the system automatically revert to manual mode (conventional pump therapy using programmed basal rates —

If the "time to exit" expires without a BG entry, the pump will revert to manual

User must enter a BG value into the pump to return to SmartGuard following

If there is no CGM data \geq 20 min, the pump will revert to manual mode.

When CGM data returns, Control-IQ will automatically turn back on.

If there is an "Automated Delivery Restriction" alarm, the user will be

mode (the Pod will not return to automated mode on its own).

Loop automation resumes when CGM data returns.

There is no option for manual mode in the iLet.

resume insulin delivery once the CGM is re-connected.

prompted to confirm CGM accuracy, and then will have to switch to manual mode. The user must switch back to automated mode after 5 min in manual

When there is no CGM data > 15 min, the pump will revert to manual mode.

After 72 hours in BG-run mode, the iLet can no longer deliver insulin. It will

The iLet will display total daily insulin dose, basal insulin and meal insulin

doses, which could be used to inform multiple daily injection doses, if needed.

adjusted to impact automated

than usual/"Less" than usual).

as needed in the bolus menu.

User completes a meal

as needed in the bolus menu / bolus

as needed in the bolus menu / bolus

as needed in the bolus menu / bolus

give?

calculator.

calculator.

calculator.

bolus menu.

Bolus automation?

Diabetes Technology. Deciphered. PANTHER program.org

MiniMed[™] 780G

Basal automation?

"Auto Basal" calculated from total daily

insulin, which is updated each day at

5 min based on recent CGM glucose

Increases or decreases the

glucose range.

programmed basal rates every 5

"Adaptive Basal" calculated from

total daily insulin, which is updated at each Pod change. Adaptive Basal

is adjusted every 5 min based on a

60 min prediction of CGM glucose,

aiming for the target glucose value.

Adjusts basal rates every 5 minutes

based on a 6-hour predicted CGM

glucose, aiming for the correction

target (middle of correction range).

Insulin Automation is initialized by

delivery adjusts every 5 minutes

entering user's weight. Basal insulin

based on CGM glucose trends and

analysis of the user's daily glucose

patterns.

Can users

(programmed

basal rates

in Smart-

Guard)

Yes

No

(programmed

basal rates

Mode)

Yes

N/A

pump settings.

correction boluses.

(manual mode).

for two reasons:

CGM data returns.

basal insulin needs.

boluses from the iLet.

on the CGM trend arrow.

Pre-bolus for all meals and snacks.

are not used in Automated

are not used

adjust

basal

rates?

No

adapts over time based on the iLet's

Can users

adjust I:C

ratios?

Yes

Yes

Yes

Yes

of CGM glucose, aiming for the target

midnight. Auto Basal is adjusted every

trends, aiming for the target glucose value. minutes based on a 30 min prediction

t:slim X2[™] & Mobi Omnipod® 5

twiist[™] AID System

powered by Tidepool

iLet Bionic Pancreas

When using AID:

MiniMed[™] 780G

t:slim X2[™] & Mobi

Omnipod® 5

twiist™ AID System

powered by Tidepool

iLet Bionic Pancreas

MiniMed[™] 780G

t:slim X2[™] & Mobi

Omnipod® 5

twiist™ AID System

powered by Tidepool

iLet Bionic Pancreas

All Pumps

MiniMed[™] 780G

t:slim X2[™] & Mobi

Omnipod® 5

twiist™ AID System

powered by Tidepool

iLet Bionic Pancreas

MiniMed[™] 780G

t:slim X2[™] & Mobi

Omnipod® 5

twiist™ AID System

powered by Tidepool

iLet Bionic Pancreas